

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	49.2	
2 - 1	53.1	
4 - 1	50.6	
5 - 1	48.2	
8 - 1	50.2	
9 - 1	52.5	
11 - 1	55.2	
12 - 1	53.4	
13 - 1	46.3	
14 - 1	49.2	
15 - 1	53.7	
16 - 1	47.5	
18 - 1	51.9	
19 - 1	48.3	

Westgard	تفسیر	IgM	آزمایش	BS-480 اتوآنالایزر	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	50.94	SD	3.92	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	50.66	SD	2.68	CV	5.29			

