

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	319.14
SD	8.66
CV	2.71
Day-Run	Value
1 - 1	316.4
2 - 1	321
5 - 1	326
6 - 1	308.6
8 - 1	330
11 - 1	315
11 - 2	315
13 - 1	314
14 - 1	310.5
15 - 1	320
16 - 1	300
18 - 1	320
19 - 1	322
21 - 1	330
22 - 1	329
23 - 1	325
25 - 1	317.2
27 - 1	336
28 - 1	317
29 - 1	310

High - CV : 1.74	
Mean	416.72
SD	7.24
CV	1.74
Day-Run	Value
1 - 1	411
2 - 1	420
5 - 1	411
6 - 1	419.4
8 - 1	406
11 - 1	436
11 - 2	420
13 - 1	411
14 - 1	418
15 - 1	404.3
16 - 1	425
18 - 1	430
19 - 1	428
20 - 1	420
21 - 1	410
23 - 1	415
25 - 1	411
27 - 1	422
28 - 1	423
29 - 1	413

Westgard	تفسیر	LDH	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Normal Control	biorex	Mean	312.94	SD	8.52	Exp. Date	1403/1/17	Lot. No	035

