

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	132.9	
2 - 1	139.4	
5 - 1	142	
6 - 1	145	1:2S
8 - 1	139	
9 - 1	141	
11 - 1	136	
12 - 1	137.4	
13 - 1	140	
14 - 1	137	
15 - 1	135.8	
16 - 1	138.3	
18 - 1	140	
19 - 1	138	
20 - 1	142	
21 - 1	141	
22 - 1	142	
23 - 1	143	4:1S
25 - 1	136.7	
27 - 1	141	
28 - 1	139	
29 - 1	139.5	

Westgard	تفسیر	LDL	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	138.2	SD	2.62	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	139.36	SD	2.76	CV	1.98			

