

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
11 - 1	11.9	
12 - 1	11.9	
13 - 1	12.1	
14 - 1	11.6	
15 - 1	11.59	
16 - 1	11.96	
18 - 1	11.4	
19 - 1	11.72	
20 - 1	11.9	
21 - 1	12.1	
22 - 1	11.3	
23 - 1	11.6	
25 - 1	11.9	
27 - 1	11.58	
28 - 1	11.6	
30 - 1	10.98	

Westgard	تفسیر	Pb	آزمایش	Atomic Absorption Spectrometry	دستگاه	9	ماه	1402	سال
Control	آزمایشگاه	Mean	11.47	SD	0.41	Exp. Date	1402/11/30	Lot. No	
Next Month	Mean	11.7	SD	0.3	CV	2.59			

