

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	8.33	
2 - 1	8.45	
4 - 1	8.38	
5 - 1	9.22	
6 - 1	9	
7 - 1	10.2	
8 - 1	9.76	
9 - 1	8.6	
11 - 1	9.38	
12 - 1	9.5	
13 - 1	10.25	
14 - 1	9.2	
15 - 1	10	
16 - 1	8.16	1:2S
18 - 1	9.86	
19 - 1	10.22	1:3S
20 - 1	9.75	
21 - 1	10.5	
22 - 1	9.5	
23 - 1	8.9	
25 - 1	10.3	
27 - 1	10	
29 - 1	9.85	
30 - 1	9.15	

Westgard	تفسیر	Progesterone	آزمایش	SIEMENS Centaur CP	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	9.52	SD	0.62	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	9.44	SD	0.7	CV	7.41			

