

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	147.2	
2 - 1	147.6	
4 - 1	148	
5 - 1	147.2	
8 - 1	147.6	
9 - 1	148.2	
11 - 1	148.5	
12 - 1	147.2	
13 - 1	147	
14 - 1	148	
15 - 1	148.2	
16 - 1	147.9	
18 - 1	148.3	
19 - 1	148	
20 - 1	147.4	
21 - 1	148.2	
22 - 1	147.9	
23 - 1	146.5	
25 - 1	147.7	

Westgard	تفسیر	Sodium (Na)	آزمایش	AUDICOM	دستگاه	9	ماه	1402	سال
Control	زیست شیمی	Mean	147.5	SD	0.53	Exp. Date	1404/1/14	Lot. No	2103061-4
Next Month	Mean	147.72	SD	0.52	CV	0.35			

