

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	2.2	
2 - 1	2.2	
4 - 1	2.42	
5 - 1	2.07	1:2S
6 - 1	2.37	
7 - 1	2.39	
8 - 1	2.67	
9 - 1	2.52	
11 - 1	2.6	
12 - 1	2.53	
14 - 1	2.85	1:2S
15 - 1	2.45	
16 - 1	2.5	
18 - 1	2.42	
19 - 1	2.27	
20 - 1	2.54	
21 - 1	2.2	
22 - 1	2.42	
23 - 1	2.35	
25 - 1	2.25	
27 - 1	2.66	
28 - 1	2.73	
29 - 1	2.41	
30 - 1	2.25	

Westgard	تفسیر	T3	آزمایش	SIEMENS Centaur CP	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	2.47	SD	0.18	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	2.43	SD	0.19	CV	7.83			

