

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	12.2	
2 - 1	11.2	
4 - 1	11	
5 - 1	12.8	1:2S
6 - 1	12	
7 - 1	12.7	
8 - 1	12.5	
9 - 1	12.4	
11 - 1	10.9	
12 - 1	11.1	
14 - 1	11.7	
15 - 1	12.5	
16 - 1	11.1	
18 - 1	11.2	
19 - 1	10.9	
20 - 1	11.95	
21 - 1	12.1	
22 - 1	12	
23 - 1	11.5	
25 - 1	11.7	
27 - 1	12.2	
28 - 1	12.6	
29 - 1	12	

Westgard	تفسیر	T4	آزمایش	SIEMENS Centaur CP	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	11.48	SD	0.65	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	11.84	SD	0.62	CV	5.21			

