

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	302.5	
2 - 1	304.7	
5 - 1	305	
6 - 1	311	
8 - 1	304	1:2S
9 - 1	301	1:2S 2:2S
11 - 1	297.3	
12 - 1	307.1	
13 - 1	303.7	
14 - 1	312.1	
15 - 1	298.1	
16 - 1	301.4	
18 - 1	302	
19 - 1	303	
20 - 1	299	
21 - 1	300	
22 - 1	303	
23 - 1	305	
25 - 1	311	1:3S
27 - 1	300	
28 - 1	306.4	
29 - 1	314	

Westgard	تفسیر	TIBC	آزمایش	اتوآنالایزر BS-800	دستگاه	9	ماه	1402	سال
Control	biorex	Mean	302.72	SD	6.1	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	304.15	SD	4.59	CV	1.51			

