

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	60.14
SD	3.39
CV	5.64
Day-Run	Value
23 - 1	55.8
24 - 1	59.4
25 - 1	58.8
26 - 1	64.8
28 - 1	61.9
29 - 1	68.2

High - CV : 4.41	
Mean	111.67
SD	4.93
CV	4.41
Day-Run	Value
23 - 1	104
24 - 1	116
25 - 1	111
26 - 1	110
28 - 1	111
29 - 1	118

Westgard	تفسیر	AFP	آزمایش	2000 ایمولایت XPI-	دستگاه	11	ماه	1402	سال
Normal Control	Randox	Mean	59.25	SD	2.86	Exp. Date	1403/11/24	Lot. No	8374 MS

