

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	4.16
SD	0.05
CV	1.18
Day-Run	Value
2 - 1	4.2
3 - 1	4.2
4 - 1	4.2
5 - 1	4.1
6 - 1	4.2
7 - 1	4.1
9 - 1	4.2
10 - 1	4.1
11 - 1	4.2
12 - 1	4.2
13 - 1	4.2
14 - 1	4.2
16 - 1	4.2
17 - 1	4.1
18 - 1	4.1
19 - 1	4.1
20 - 1	4.1
21 - 1	4.1
23 - 1	4.2
24 - 1	4.2
25 - 1	4.2
26 - 1	4.1
27 - 1	4.2
28 - 1	4.2
30 - 1	4.2

High - CV : 1.38	
Mean	5.11
SD	0.07
CV	1.38
Day-Run	Value
2 - 1	5.1
3 - 1	5.2
4 - 1	5.2
5 - 1	5
6 - 1	5.2
7 - 1	5.1
9 - 1	5.1
10 - 1	5.2
11 - 1	5.1
12 - 1	5
13 - 1	5.1
14 - 1	5
16 - 1	5
17 - 1	5.1
18 - 1	5.1
19 - 1	5.1
20 - 1	5.1
21 - 1	5.2
23 - 1	5.1
24 - 1	5
25 - 1	5.1
26 - 1	5.1
27 - 1	5.2
28 - 1	5.2
30 - 1	5.1

Westgard	تفسیر	Albomin	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	4.17	SD	0.08	Exp. Date	1403/1/17	Lot. No	035

