

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.24	
Mean	29.14	Mean	138.38
SD	1.25	SD	3.1
CV	4.28	CV	2.24
Day-Run	Value	Day-Run	Value
1 - 1	29.2	1 - 1	144
2 - 1	28.5	2 - 1	139.6
3 - 1	27.1	3 - 1	143.5
4 - 1	27.6	4 - 1	141.2
7 - 1	31	7 - 1	136.7
8 - 1	31	8 - 1	136
9 - 1	28	9 - 1	138
10 - 1	28	10 - 1	134.9
11 - 1	30.7	11 - 1	141
12 - 1	30	12 - 1	133.1
14 - 1	28.8	14 - 1	136.4
15 - 1	29.2	15 - 1	135
16 - 1	30	16 - 1	137
17 - 1	29.1	17 - 1	136.6
18 - 1	29.3	18 - 1	140
21 - 1	28.8	21 - 1	137
23 - 1	28.4	23 - 1	138
24 - 1	28	24 - 1	136
25 - 1	30.7	25 - 1	137.2
28 - 1	27.5	28 - 1	143.2
29 - 1	31	29 - 1	136.9
30 - 1	31	30 - 1	143

Westgard	تفسیر	SGPT (ALT)	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	28.6	SD	1.5	Exp. Date	1403/1/17	Lot. No	035

