

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	28.6
SD	0.94
CV	3.29
Day-Run	Value
2 - 1	30.1
3 - 1	28.7
4 - 1	28.9
5 - 1	29.1
6 - 1	29.8
7 - 1	29.8
9 - 1	28.4
10 - 1	30
11 - 1	28.4
12 - 1	28.2
13 - 1	27.8
14 - 1	27.3
16 - 1	27
17 - 1	29.4
18 - 1	27.4
19 - 1	29.5
20 - 1	27.3
21 - 1	27.2
23 - 1	28
24 - 1	29.6
25 - 1	28.6
26 - 1	28
27 - 1	29.2
28 - 1	28.6
30 - 1	28.7

High - CV : 1.79	
Mean	137.38
SD	2.46
CV	1.79
Day-Run	Value
2 - 1	134.5
3 - 1	134
4 - 1	138
5 - 1	133
6 - 1	138
7 - 1	136
9 - 1	134.8
10 - 1	135.4
11 - 1	132.8
12 - 1	137.3
13 - 1	138
14 - 1	139
16 - 1	140
17 - 1	137
18 - 1	138
19 - 1	139
20 - 1	138
21 - 1	135
23 - 1	138
24 - 1	141.3
25 - 1	142.1
26 - 1	140.1
27 - 1	137.1
28 - 1	138
30 - 1	140

Westgard	تفسیر	SGPT ( ALT )	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	29	SD	1.5	Exp. Date	1403/1/17	Lot. No	035

