

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	38.39
SD	0.77
CV	1.99
Day-Run	Value
2 - 1	39.2
3 - 1	37
4 - 1	38
5 - 1	38.2
6 - 1	39
7 - 1	40
9 - 1	39
10 - 1	37.8
11 - 1	38.8
12 - 1	37.7
13 - 1	38.1
14 - 1	38.2
16 - 1	39.2
17 - 1	38
18 - 1	38.1
19 - 1	38.2
20 - 1	36.9
21 - 1	37.8
23 - 1	38
24 - 1	39.4
25 - 1	39.1
26 - 1	39.1
27 - 1	37.5
28 - 1	39
30 - 1	38.4

High - CV : 1.37	
Mean	146.04
SD	2
CV	1.37
Day-Run	Value
2 - 1	145
3 - 1	146.3
4 - 1	142
5 - 1	144
6 - 1	146
7 - 1	146.3
9 - 1	146
10 - 1	143
11 - 1	145
12 - 1	146
13 - 1	147
14 - 1	149
16 - 1	146
17 - 1	148
18 - 1	147
19 - 1	145
20 - 1	150
21 - 1	147
23 - 1	144
24 - 1	148
25 - 1	144.2
26 - 1	148
27 - 1	143.2
28 - 1	146
30 - 1	149

Westgard	تفسیر	SGOT ( AST )	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	38	SD	1	Exp. Date	1403/1/17	Lot. No	035

