

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	79.6	
2 - 1	79.5	
3 - 1	77.4	
4 - 1	79	
7 - 1	77.3	
8 - 1	79.5	
9 - 1	80.3	
10 - 1	82.7	
11 - 1	80	
12 - 1	83.4	
14 - 1	80	
15 - 1	80	
16 - 1	79	
17 - 1	80.2	
18 - 1	83	
21 - 1	81	
23 - 1	82.7	
24 - 1	80	
25 - 1	79	
28 - 1	81	
29 - 1	75.8	1:2S
30 - 1	82.3	1:2S 2:2S

Westgard	تفسیر	Amylase	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	80.76	SD	1.94	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	79.84	SD	2.05	CV	2.56			

