

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	81	
2 - 1	81	
3 - 1	79.2	
4 - 1	82.9	
5 - 1	81.8	
6 - 1	83	
7 - 1	78.5	
9 - 1	83.5	
10 - 1	83	
11 - 1	81	
12 - 1	79.7	
13 - 1	81.4	
14 - 1	79.9	
16 - 1	81.9	
17 - 1	80.1	
18 - 1	82.3	
18 - 1	82.3	
19 - 1	79.3	
20 - 1	81.5	
21 - 1	83	
23 - 1	81.2	
24 - 1	80.2	
25 - 1	81.2	
26 - 1	80.3	
27 - 1	82	
28 - 1	77.1	
28 - 1	77.1	
30 - 1	75.8	1:2S

Westgard	تفسیر	Amylase	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	80.33	SD	1.66	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	80.76	SD	1.94	CV	2.4			

