

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	24.3	
3 - 1	25.3	
4 - 1	26.1	
5 - 1	24.8	
6 - 1	23.8	
7 - 1	24.6	
9 - 1	25.7	
10 - 1	21.1	1:2S
11 - 1	22.5	
12 - 1	23.1	
13 - 1	24.5	
14 - 1	23.9	
16 - 1	25.1	
17 - 1	25.3	
18 - 1	23.4	
19 - 1	24.9	
20 - 1	24.1	
21 - 1	25.7	
23 - 1	23.2	
24 - 1	23.9	
25 - 1	23.4	
26 - 1	24.8	
27 - 1	22.9	
28 - 1	23.7	
30 - 1	25.1	

Westgard	تفسیر	BHCG	آزمایش	2000 ایمولایت XPI-	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	24.42	SD	1.29	Exp. Date	1403/8/30	Lot. No	
Next Month	Mean	24.21	SD	1.15	CV	4.76			

