

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	17.4	
2 - 1	18	
3 - 1	17.9	
4 - 1	17.7	
7 - 1	18.1	
8 - 1	18	
9 - 1	17.9	
10 - 1	17.63	
12 - 1	18.2	
14 - 1	17.5	
15 - 1	17.6	
16 - 1	18.1	
17 - 1	18.1	
18 - 1	17.8	1:3S
21 - 1	18	
23 - 1	17.9	
24 - 1	18.3	
25 - 1	18	
28 - 1	17.6	
29 - 1	17.7	
30 - 1	17.9	

Westgard	تفسیر	C4	آزمایش	BS-480 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	17.87	SD	0.27	Exp. Date	1402/11/30	Lot. No	641-01-030/2
Next Month	Mean	17.87	SD	0.24	CV	1.34			

