

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	48.5	
3 - 1	48.1	
4 - 1	52.8	
5 - 1	52.6	
6 - 1	50.7	
7 - 1	51.9	
9 - 1	52.3	
10 - 1	48.8	
11 - 1	48.9	
12 - 1	49.1	
13 - 1	50.8	
14 - 1	51.9	
16 - 1	52.1	
17 - 1	51.4	
18 - 1	52.7	
19 - 1	52.9	
20 - 1	52.5	
21 - 1	53.1	
23 - 1	53.8	
24 - 1	50.4	
25 - 1	50.9	
26 - 1	52.1	
27 - 1	53.4	
28 - 1	52.9	
30 - 1	51.2	

Westgard	تفسیر	CA 125	آزمایش	2000 ایمولایت XPI-	دستگاه	10	ماه	1402	سال
Control	Seronorm	Mean	52.03	SD	2.66	Exp. Date	1403/8/30	Lot. No	2011909
Next Month	Mean	51.43	SD	1.65	CV	3.21			

