

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

| Day-Run | Value | Description |
|---------|-------|-------------|
| 1 - 1 | 52.2 | |
| 2 - 1 | 52.9 | |
| 3 - 1 | 46 | 1:2S |
| 4 - 1 | 53.8 | |
| 7 - 1 | 46.3 | 1:2S |
| 14 - 1 | 48.5 | |
| 15 - 1 | 49.5 | |
| 16 - 1 | 57.6 | |
| 17 - 1 | 52.1 | |
| 18 - 1 | 57.1 | |
| 23 - 1 | 45.8 | 1:2S |
| 25 - 1 | 52.1 | |
| 26 - 1 | 48.2 | |
| 28 - 1 | 53.7 | |
| 29 - 1 | 51.4 | |
| 30 - 1 | 50.7 | |

| | | | | | | | | | |
|------------|----------|---------|--------|--------------------|--------|-----------|------------|---------|---------|
| Westgard | تفسیر | CA 15-3 | آزمایش | 2000 ایمولایت XPI- | دستگاه | 11 | ماه | 1402 | سال |
| Control | Seronorm | Mean | 53.3 | SD | 2.74 | Exp. Date | 1402/11/25 | Lot. No | 2011909 |
| Next Month | Mean | 51.12 | SD | 3.6 | CV | 7.03 | | | |

