

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	22.4	
2 - 1	23.4	
3 - 1	21.9	
4 - 1	23.5	
7 - 1	23.8	
8 - 1	24.1	
9 - 1	22.9	
10 - 1	22.7	
11 - 1	23.5	
12 - 1	23.7	
14 - 1	22.2	
15 - 1	23.1	
16 - 1	23.5	
17 - 1	23.4	
18 - 1	22.8	
21 - 1	23.7	
23 - 1	22.6	
24 - 1	23.1	
25 - 1	23.8	
25 - 2	22.3	
26 - 1	23.4	
28 - 1	21.3	1:2S
29 - 1	22.6	
30 - 1	23.9	

Westgard	تفسیر	CEA	آزمایش	2000 ایمولایت XPI-	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	23.45	SD	0.83	Exp. Date	1402/11/15	Lot. No	36502028
Next Month	Mean	23.07	SD	0.71	CV	3.06			

