

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 0.75	
Mean	173.1	Mean	264.43
SD	2.47	SD	1.99
CV	1.43	CV	0.75
Day-Run	Value	Day-Run	Value
1 - 1	177	1 - 1	267
2 - 1	177	2 - 1	262
3 - 1	172	3 - 1	262
4 - 1	176	4 - 1	263
7 - 1	172	7 - 1	265
8 - 1	170	8 - 1	268
9 - 1	174	9 - 1	264
10 - 1	170	10 - 1	265
11 - 1	175	11 - 1	263
12 - 1	170	12 - 1	267
14 - 1	172	14 - 1	266
15 - 1	177	15 - 1	260
16 - 1	176	16 - 1	264
18 - 1	172	18 - 1	264
21 - 1	172	21 - 1	264
23 - 1	173	23 - 1	263
24 - 1	172	24 - 1	265
25 - 1	171	25 - 1	266
28 - 1	170	28 - 1	265
29 - 1	175	29 - 1	263
30 - 1	172	30 - 1	267

Westgard	تفسیر	Cholesterol	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	175.24	SD	2.24	Exp. Date	1403/1/17	Lot. No	035

