

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	48	
3 - 1	48.2	
4 - 1	48.7	1:2S
5 - 1	46.8	
6 - 1	47.5	
7 - 1	48	
9 - 1	46.9	
10 - 1	48	
11 - 1	47	
12 - 1	47.3	
13 - 1	47.6	
14 - 1	46.5	
16 - 1	47.7	
17 - 1	47.2	
18 - 1	48.5	
19 - 1	48	
20 - 1	48.5	
21 - 1	47.3	
23 - 1	48.5	
24 - 1	47.2	
25 - 1	48.2	
26 - 1	47.3	
27 - 1	46.2	
28 - 1	47.3	
30 - 1	47.8	

Westgard	تفسیر	CRP-quantitative	آزمایش	BS-480 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	47.42	SD	0.62	Exp. Date	1403/6/8	Lot. No	326-2-034/3
Next Month	Mean	47.61	SD	0.66	CV	1.38			

