

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	24.9	
2 - 1	25.1	
3 - 1	25.3	
4 - 1	24.6	
7 - 1	24.9	
8 - 1	25.2	
9 - 1	25.1	
10 - 1	24.6	
12 - 1	24.5	
15 - 1	24.6	
21 - 1	24.8	
24 - 1	25.1	
28 - 1	25	
30 - 1	25.3	

Westgard	تفسیر	Cu,urine	آزمایش	Atomic Absorption Spectrometry	دستگاه	11	ماه	1402	سال
Control	آزمایشگاه	Mean	24.99	SD	0.29	Exp. Date	1403/1/30	Lot. No	
Next Month	Mean	24.93	SD	0.27	CV	1.09			

