

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	30.4	
2 - 1	32	
3 - 1	29.9	
4 - 1	31.1	
7 - 1	31.2	
8 - 1	32	
9 - 1	30.3	
10 - 1	32.1	
12 - 1	29.7	
14 - 1	30.1	
15 - 1	31.2	
16 - 1	31.8	
17 - 1	31	
18 - 1	29.4	
21 - 1	32	
23 - 1	32	
24 - 1	31.6	
25 - 1	30.1	
28 - 1	32	
29 - 1	32	
30 - 1	30.4	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	30.68	SD	0.75	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	31.06	SD	0.91	CV	2.92			

