

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	29.8	
3 - 1	31.6	
4 - 1	30.8	
5 - 1	29.5	
6 - 1	29.6	1:2S
7 - 1	31.7	
9 - 1	30.9	
10 - 1	31	
11 - 1	31.3	
12 - 1	30.1	
13 - 1	30.6	
14 - 1	31.5	
16 - 1	31.4	
18 - 1	31.7	4:1S
17 - 1	30.2	
19 - 1	31.8	
20 - 1	29.9	
21 - 1	30.4	
23 - 1	29.8	
24 - 1	30.1	
25 - 1	31.2	
26 - 1	31.3	
27 - 1	29.7	
28 - 1	30.2	
30 - 1	31	

Westgard	تفسیر	Ceruloplasmin	آزمایش	BS-480 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	30.37	SD	0.76	Exp. Date	1406/9/14	Lot. No	641-01-030/2
Next Month	Mean	30.68	SD	0.75	CV	2.45			

