

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	128	
3 - 1	128	
4 - 1	128	
5 - 1	128	
6 - 1	126	
7 - 1	126	
9 - 1	128	
9 - 1	128	
10 - 1	127	
11 - 1	127	
12 - 1	129	
13 - 1	127	
14 - 1	127	
16 - 1	127	
17 - 1	128	
18 - 1	127	
19 - 1	129	
20 - 1	128	
21 - 1	130	1:2S
23 - 1	128.2	
24 - 1	126.5	
25 - 1	128	
26 - 1	129	
27 - 1	126.6	
28 - 1	128	
30 - 1	128	

Westgard	تفسیر	Copper ( Cu )	آزمایش	اتوانالایزر BS-800	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	127.58	SD	1.08	Exp. Date	1402/11/30	Lot. No	hsn.35
Next Month	Mean	127.7	SD	0.96	CV	0.75			

