

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	2.65
SD	0.04
CV	1.41
Day-Run	Value
2 - 1	2.66
3 - 1	2.65
4 - 1	2.64
5 - 1	2.63
6 - 1	2.6
7 - 1	2.65
9 - 1	2.7
10 - 1	2.7
11 - 1	2.7
12 - 1	2.6
13 - 1	2.7
14 - 1	2.64
16 - 1	2.66
17 - 1	2.7
18 - 1	2.64
19 - 1	2.67
20 - 1	2.6
21 - 1	2.66
23 - 1	2.7
24 - 1	2.6
25 - 1	2.63
26 - 1	2.61
27 - 1	2.6
28 - 1	2.6
30 - 1	2.64

High - CV : 1.71	
Mean	3.73
SD	0.06
CV	1.71
Day-Run	Value
2 - 1	3.73
3 - 1	3.76
4 - 1	3.78
5 - 1	3.7
6 - 1	3.8
7 - 1	3.7
9 - 1	3.9
10 - 1	3.8
11 - 1	3.8
12 - 1	3.67
13 - 1	3.7
14 - 1	3.65
16 - 1	3.77
17 - 1	3.7
18 - 1	3.71
19 - 1	3.8
20 - 1	3.68
21 - 1	3.7
23 - 1	3.7
24 - 1	3.6
25 - 1	3.7
26 - 1	3.8
27 - 1	3.69
28 - 1	3.7
30 - 1	3.76

Westgard	تفسیر	Creatinine	آزمایش	BS-800 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	2.65	SD	0.04	Exp. Date	1403/1/17	Lot. No	035

