

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	0.47
SD	0.03
CV	6.88
Day-Run	Value
1 - 1	0.43
2 - 1	0.55
3 - 1	0.46
4 - 1	0.48
7 - 1	0.45
8 - 1	0.46
9 - 1	0.48
10 - 1	0.46
11 - 1	0.52
12 - 1	0.53
14 - 1	0.47
15 - 1	0.5
16 - 1	0.43
17 - 1	0.47
17 - 1	0.47
18 - 1	0.47
21 - 1	0.46
23 - 1	0.47
24 - 1	0.41
25 - 1	0.43
28 - 1	0.47
29 - 1	0.46
30 - 1	0.48

High - CV : 4.37	
Mean	2.22
SD	0.1
CV	4.37
Day-Run	Value
1 - 1	2.05
2 - 1	2.36
3 - 1	2.3
4 - 1	2.27
7 - 1	2.05
8 - 1	2.3
9 - 1	2.3
10 - 1	2.3
11 - 1	2.25
12 - 1	2.15
14 - 1	2.32
15 - 1	2.2
16 - 1	2.3
17 - 1	2.18
18 - 1	2.3
21 - 1	2.1
23 - 1	2.2
24 - 1	2.3
25 - 1	2.1
28 - 1	2.2
29 - 1	2.3
30 - 1	2.08

Westgard	تفسیر	D.Bili	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	0.47	SD	0.03	Exp. Date	1403/1/17	Lot. No	035

