

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	73.59
SD	1.18
CV	1.61
Day-Run	Value
1 - 1	76
2 - 1	74
3 - 1	71
4 - 1	72
7 - 1	72
8 - 1	74
9 - 1	74
10 - 1	74
11 - 1	73
12 - 1	73
14 - 1	75
15 - 1	73
16 - 1	73
17 - 1	74
18 - 1	72
21 - 1	74
23 - 1	74
24 - 1	74
25 - 1	75
28 - 1	75
29 - 1	74
30 - 1	73

High - CV : 1.29	
Mean	229.32
SD	2.97
CV	1.29
Day-Run	Value
1 - 1	234
2 - 1	232
3 - 1	229
4 - 1	231
7 - 1	230
8 - 1	220
9 - 1	230
10 - 1	230
11 - 1	226
12 - 1	233
14 - 1	228
15 - 1	225
16 - 1	230
17 - 1	230
18 - 1	228
21 - 1	229
23 - 1	233
24 - 1	230
25 - 1	229
28 - 1	230
29 - 1	228
30 - 1	230

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	73.6	SD	0.91	Exp. Date	1403/1/17	Lot. No	035

