

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High - CV : 1.55	
Mean	73.6	Mean	229.2
SD	0.91	SD	3.56
CV	1.24	CV	1.55
Day-Run	Value	Day-Run	Value
2 - 1	74	2 - 1	229
3 - 1	74	3 - 1	232
4 - 1	75	4 - 1	233
5 - 1	73	5 - 1	230
6 - 1	74	6 - 1	232
7 - 1	74	7 - 1	228
9 - 1	74	9 - 1	227
10 - 1	75	10 - 1	221
11 - 1	73	11 - 1	230
12 - 1	73	12 - 1	230
13 - 1	75	13 - 1	220
14 - 1	72	14 - 1	225
16 - 1	74	16 - 1	231
17 - 1	72	17 - 1	225
18 - 1	73	18 - 1	232
19 - 1	74	19 - 1	231
20 - 1	73	20 - 1	227
21 - 1	74	21 - 1	226
23 - 1	74	23 - 1	233
24 - 1	75	24 - 1	231
25 - 1	74	25 - 1	233
26 - 1	73	26 - 1	232
27 - 1	73	27 - 1	230
28 - 1	73	28 - 1	230
30 - 1	72	30 - 1	232

Westgard	تفسیر	FBS	آزمایش	BS-800 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	73.33	SD	1.08	Exp. Date	1403/1/17	Lot. No	035

