

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	23.1	
3 - 1	24.8	
4 - 1	24.1	
5 - 1	24.9	
6 - 1	23.8	
7 - 1	25.6	
9 - 1	26.1	1:2S
10 - 2	23.5	
11 - 1	23.9	
12 - 1	23.8	
13 - 1	24.1	
14 - 1	26.2	1:2S
16 - 1	25.2	
17 - 1	25.7	
18 - 1	24.9	
19 - 1	23.7	
20 - 1	24.1	
21 - 1	24.5	
23 - 1	24.8	
24 - 1	24.1	
25 - 1	23.9	
26 - 1	24.6	
27 - 1	24.8	
28 - 1	23.7	
30 - 1	25.1	

Westgard	تفسیر	FSH	آزمایش	2000 ایمولایت XPI-	دستگاه	10	ماه	1402	سال
Control	آزمایشگاه	Mean	24.18	SD	0.8	Exp. Date	1403/8/30	Lot. No	
Next Month	Mean	24.52	SD	0.82	CV	3.34			

