

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	155.7	
3 - 1	148.2	
4 - 1	144.6	
5 - 1	142.9	
6 - 1	147.8	
7 - 1	153.2	
9 - 1	150.7	
10 - 1	153.7	
11 - 1	149.8	
12 - 1	141.5	
13 - 1	147.6	
14 - 1	140	1:3S Cusum Reject
16 - 1	151	
17 - 1	152.7	
18 - 1	151.5	
19 - 1	149.2	
20 - 1	149.8	
21 - 1	152.1	
23 - 1	155.6	
24 - 1	157.6	
25 - 1	139.7	
26 - 1	142.7	
27 - 1	154.2	
28 - 1	136.5	
30 - 1	139.2	

Westgard	تفسیر	Ferritin	آزمایش	SIEMENS Centaur CP	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	144.84	SD	6.56	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	148.3	SD	5.88	CV	3.97			

