

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	7.35	
2 - 1	7.3	
3 - 1	7.68	
4 - 1	7.35	
7 - 1	7.22	
8 - 1	7.39	
9 - 1	7.94	1:2S
10 - 1	7.61	
11 - 1	7.67	
12 - 1	7.75	
14 - 1	7.25	
15 - 1	7.4	
16 - 1	7.1	
17 - 1	7.35	
18 - 1	7.3	
21 - 1	7.6	
23 - 1	7.8	
24 - 1	7.63	
25 - 1	7.42	
26 - 1	7.24	
28 - 1	7.44	
29 - 1	7.1	
30 - 1	7.32	

Westgard	تفسیر	Free T3	آزمایش	SIEMENS Centaur CP	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	7.42	SD	0.2	Exp. Date	1402/11/25	Lot. No	36-02-026
Next Month	Mean	7.44	SD	0.23	CV	3.02			

