

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	7.25	
3 - 1	7.12	
4 - 1	7.27	
5 - 1	7.53	
6 - 1	7.65	
7 - 1	7.15	
9 - 1	7.24	
10 - 1	7.6	
11 - 1	7.55	
12 - 1	7.73	
13 - 1	7.43	
14 - 1	7.52	
16 - 1	7.22	
17 - 1	7	1:2S
18 - 1	7.25	
19 - 1	7.22	
20 - 1	7.36	
21 - 1	7.75	
23 - 1	7.5	
24 - 1	7.65	
25 - 1	7.46	
26 - 1	7.35	
27 - 1	7.66	
28 - 1	7.45	
30 - 1	7.5	

Westgard	تفسیر	Free T3	آزمایش	SIEMENS Centaur CP	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	7.42	SD	0.17	Exp. Date	1403/8/30	Lot. No	36-02-026
Next Month	Mean	7.42	SD	0.2	CV	2.75			

