

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	70.9	
3 - 1	68	
4 - 1	70.1	
5 - 1	66.1	
6 - 1	69	
7 - 1	68	
9 - 1	68.3	
10 - 1	68.1	
11 - 1	69.7	
12 - 1	67.9	
13 - 1	65.3	
14 - 1	68.2	
16 - 1	67.6	
17 - 1	67.1	
18 - 1	68.2	
19 - 1	65.1	
20 - 1	65.6	
21 - 1	65.4	
23 - 1	69.1	
24 - 1	66.1	
25 - 1	67.1	
26 - 1	68.2	
27 - 1	66.1	
28 - 1	63.4	1:2S
30 - 1	68.4	

Westgard	تفسیر	Gamma GT	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	67.71	SD	1.84	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67.48	SD	1.75	CV	2.59			

