

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	13.3	Mean	16.43	Mean	5.19
SD	0.17	SD	0.1	SD	0.04
CV	1.3	CV	0.63	CV	0.75
Day-Run	Value	Day-Run	Value	Day-Run	Value
2 - 1	13.3	2 - 1	16.4	2 - 1	5.2
3 - 1	13.2	3 - 1	16.3	3 - 1	5.1
4 - 1	13.1	4 - 1	16.4	4 - 1	5.2
5 - 1	13.1	5 - 1	16.3	5 - 1	5.2
6 - 1	13.6	9 - 1	16.6	6 - 1	5.2
9 - 1	13.4	11 - 1	16.5	9 - 1	5.2
11 - 1	13.3	12 - 1	16.5	11 - 1	5.2
12 - 1	13.2	19 - 1	16.4	12 - 1	5.25
19 - 1	13.5			19 - 1	5.2

Westgard	تفسیر	HB	آزمایش	NIHON KOHDEN سل کانترا	دستگاه	10	ماه	1402	سال
Normal Control	R&D	Mean	13.2	SD	0.21	Exp. Date	1402/9/21	Lot. No	
High Control	R&D	Mean	16.5	SD	0.26	Exp. Date	1402/9/21	Lot. No	
Low Control	R&D	Mean	5.2	SD	0.04	Exp. Date	1402/9/21	Lot. No	

