

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	70	
2 - 1	70	
3 - 1	65.3	
4 - 1	68.1	
7 - 1	67	
8 - 1	71	1:2S
9 - 1	67	
10 - 1	68	
11 - 1	70	
12 - 1	69	
14 - 1	68.8	
15 - 1	67.8	
16 - 1	67.2	
17 - 1	66.1	
18 - 1	67.5	
21 - 1	69.8	
23 - 1	68.1	
24 - 1	67.8	
25 - 1	68.3	
28 - 1	69.7	
29 - 1	67	
30 - 1	67.8	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	67.78	SD	1.37	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	68.24	SD	1.42	CV	2.09			

