

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	68.5	
3 - 1	69	
4 - 1	68.5	
5 - 1	67	
6 - 1	69	
7 - 1	65.5	
9 - 1	68.6	
10 - 1	66.8	
11 - 1	69.7	
12 - 1	67.6	
13 - 1	70.2	1:2S
14 - 1	67.1	
16 - 1	67.4	
17 - 1	68.1	
18 - 1	69.3	
19 - 1	69.7	
20 - 1	67.5	
21 - 1	68.2	
23 - 1	66.3	
24 - 1	66	
25 - 1	66.9	
26 - 1	68	
27 - 1	66.2	
28 - 1	65	1:2S
30 - 1	68.4	

Westgard	تفسیر	HDL	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	67.7	SD	1.14	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67.78	SD	1.37	CV	2.02			

