

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	140.2	
3 - 1	141.2	
4 - 1	138.9	
5 - 1	141.3	
6 - 1	142	
7 - 1	138.7	
9 - 1	145.1	
10 - 1	144.5	
11 - 1	144.7	
12 - 1	138.3	
13 - 1	141.4	
14 - 1	138.1	
16 - 1	140.1	1:3S
17 - 1	143	
18 - 1	144.5	
19 - 1	145.2	1:3S
20 - 1	138.8	
21 - 1	144.1	
23 - 1	138.5	
24 - 1	144	
25 - 1	144.3	
26 - 1	136.6	
27 - 1	140	
28 - 1	142.1	
30 - 1	140.7	

Westgard	تفسیر	Ig A	آزمایش	BS-480 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	139.84	SD	2.91	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	141.45	SD	2.6	CV	1.84			

