

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	769.3	
3 - 1	770.3	
4 - 1	775.3	
5 - 1	767.8	
6 - 1	772.1	
7 - 1	783	
9 - 1	765.3	
10 - 1	770	
11 - 1	778.3	
12 - 1	767.2	
13 - 1	772	
14 - 1	777	
16 - 1	775.4	
17 - 1	769.2	
18 - 1	767.6	
19 - 1	768	
20 - 1	773	1:2S
21 - 1	772.9	
23 - 1	761.2	
24 - 1	773.1	
25 - 1	766.1	
26 - 1	770.3	
27 - 1	773	
28 - 1	766.2	
30 - 1	773	

Westgard	تفسیر	Ig G	آزمایش	BS-480 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	768.03	SD	8.56	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	771.06	SD	4.66	CV	0.6			

