

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	52.8	
2 - 1	52.5	
3 - 1	48.6	
4 - 1	56	
7 - 1	47.3	
8 - 1	52.3	
9 - 1	54.2	
10 - 1	56.2	
12 - 1	50.4	
14 - 1	52.7	
15 - 1	49.2	
16 - 1	54.4	
17 - 1	53.1	
18 - 1	51.6	
21 - 1	47.2	
23 - 1	55	
24 - 1	56	
25 - 1	49.3	
28 - 1	48.8	
29 - 1	54.5	
30 - 1	49.4	

Westgard	تفسیر	Ig M	آزمایش	BS-480 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	51.77	SD	2.58	Exp. Date	1402/11/7	Lot. No	641-01-030/2
Next Month	Mean	51.98	SD	2.91	CV	5.6			

