

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	229.6	
2 - 1	225	
3 - 1	223.7	
4 - 1	225.9	
7 - 1	230.8	
8 - 1	237	1:2S
9 - 1	231	
10 - 1	228	
11 - 1	230	
12 - 1	230	
14 - 1	234	
15 - 1	230	
16 - 1	220	
17 - 1	231	
18 - 1	225	
21 - 1	230	
23 - 1	224	
24 - 1	226	
25 - 1	230	
28 - 1	230	
29 - 1	220	
30 - 1	225	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	227.33	SD	4.72	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	228	SD	4.18	CV	1.83			

