

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	230	
3 - 1	230	
4 - 1	223	
5 - 1	218	
6 - 1	224.8	
7 - 1	230	
9 - 1	227	
10 - 1	228.2	
11 - 1	232.4	
12 - 1	232.5	
13 - 1	233.8	
14 - 1	219.7	4:1S
16 - 1	233.5	
17 - 1	231.5	
18 - 1	234	
19 - 1	220	
20 - 1	230	
21 - 1	220.8	
23 - 1	227	
24 - 1	226	
25 - 1	227	
27 - 1	223	
28 - 1	227.6	
30 - 1	226	

Westgard	تفسیر	Iron	آزمایش	BS-800 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	225	SD	6.25	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	227.33	SD	4.72	CV	2.08			

