

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	313.71
SD	10.64
CV	3.39
Day-Run	Value
1 - 1	295
2 - 1	317
3 - 1	305
4 - 1	318
7 - 1	328
8 - 1	342
9 - 1	320
10 - 1	320
12 - 1	306
14 - 1	296
15 - 1	314
16 - 1	310
17 - 1	305
18 - 1	313
21 - 1	320
23 - 1	302
24 - 1	318
25 - 1	320
28 - 1	312
29 - 1	310
30 - 1	317

High - CV : 2.11	
Mean	427.15
SD	9.03
CV	2.11
Day-Run	Value
1 - 1	417
2 - 1	427
3 - 1	425
4 - 1	430
7 - 1	420
9 - 1	414
10 - 1	421
11 - 1	443
12 - 1	425
15 - 1	440
16 - 1	430
17 - 1	436
18 - 1	432
21 - 1	431
23 - 1	436
24 - 1	440
25 - 1	420
28 - 1	422
29 - 1	424
30 - 1	410

Westgard	تفسیر	LDH	آزمایش	اتوانالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	317.26	SD	8.26	Exp. Date	1403/1/17	Lot. No	035

