

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	317.26
SD	8.26
CV	2.6
Day-Run	Value
2 - 1	314
3 - 1	319
4 - 1	330
5 - 1	304
6 - 1	328
7 - 1	305
9 - 1	313
10 - 1	318
11 - 1	324
12 - 1	318
13 - 1	305
14 - 1	312
16 - 1	316
17 - 1	330
18 - 1	314
19 - 1	318
20 - 1	330
21 - 1	327
23 - 1	320
24 - 1	303.2
25 - 1	321.1
26 - 1	323
27 - 1	318.2
28 - 1	311
30 - 1	310

High - CV : 1.89	
Mean	422.58
SD	7.99
CV	1.89
Day-Run	Value
2 - 1	420
3 - 1	431
4 - 1	415
5 - 1	420
6 - 1	430
7 - 1	423
9 - 1	424
10 - 1	428
11 - 1	413
12 - 1	430
13 - 1	424
14 - 1	435
16 - 1	411
17 - 1	424
18 - 1	437
19 - 1	420
20 - 1	416
21 - 1	438
23 - 1	429.3
24 - 1	415
25 - 1	417.3
26 - 1	420
27 - 1	413
28 - 1	420
30 - 1	411

Westgard	تفسیر	LDH	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	319.14	SD	8.66	Exp. Date	1403/1/17	Lot. No	035

