

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	134.4	
2 - 1	137	
3 - 1	141	
4 - 1	135.9	
7 - 1	140	
8 - 1	130	1:2S
9 - 1	137.7	
10 - 1	131	1:2S
11 - 1	133	
12 - 1	136	
14 - 1	142	
15 - 1	140	
16 - 1	138	
17 - 1	138	
21 - 1	132	
23 - 1	136	
24 - 1	130	1:2S
25 - 1	135.4	
28 - 1	137.9	
29 - 1	141	
30 - 1	138.8	

Westgard	تفسیر	LDL	آزمایش	BS-800 اتوآنالایزر	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	137.14	SD	2.88	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	136.5	SD	3.55	CV	2.6			

