

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	34.5	
2 - 1	32.1	
3 - 1	35.2	
4 - 1	34.9	
7 - 1	34.9	
8 - 1	33.5	
9 - 1	32	
10 - 1	34	
11 - 1	32.6	
12 - 1	35.4	
14 - 1	33.3	
15 - 1	34.7	
17 - 1	34.2	
18 - 1	32.1	
23 - 1	34	
24 - 1	31.8	
25 - 1	34	
26 - 1	30.9	
28 - 1	34	
29 - 1	29.5	1:2S
30 - 1	30.6	1:2S 2:2S

Westgard	تفسیر	LH	آزمایش	2000 ایمولایت XPI-	دستگاه	11	ماه	1402	سال
Control	Seronorm	Mean	33	SD	1.23	Exp. Date	1403/5/31	Lot. No	2011909
Next Month	Mean	33.2	SD	1.73	CV	5.2			

