

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	68.5	
2 - 1	69	
3 - 1	66.7	
4 - 1	65.2	
7 - 1	65	
8 - 1	68	
9 - 1	66.3	
10 - 1	66.4	
11 - 1	66.9	
12 - 1	68.8	
14 - 1	67.1	
15 - 1	67.6	
16 - 1	70	1:2S
17 - 1	67	
18 - 1	65.9	
21 - 1	65.2	
23 - 1	66.1	
24 - 1	68	
25 - 1	66	
28 - 1	65.8	
29 - 1	69.2	
30 - 1	65.2	

Westgard	تفسیر	Lipase	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	66.75	SD	1.55	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	67	SD	1.46	CV	2.18			

