

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	30.63	Mean	33.84	Mean	25.44
SD	0.37	SD	0.46	SD	0.35
CV	1.21	CV	1.35	CV	1.38
Day-Run	Value	Day-Run	Value	Day-Run	Value
2 - 1	30.6	2 - 1	34	2 - 1	25.8
3 - 1	30.8	3 - 1	34.6	3 - 1	25.2
4 - 1	30.6	4 - 1	34.2	4 - 1	25.7
5 - 1	30.5	5 - 1	34	5 - 1	25.5
6 - 1	31.3	6 - 1	33.5	6 - 1	25.3
9 - 1	31	9 - 1	34	9 - 1	25.5
11 - 1	30.4	11 - 1	33.7	11 - 1	25
12 - 1	30	12 - 1	33.6	12 - 1	26
19 - 1	30.5	19 - 1	33	19 - 1	25

Westgard	تفسیر	MCH	آزمایش	NIHON KOHDEN سل کانترا	دستگاه	10	ماه	1402	سال
Normal Control	R&D	Mean	30.46	SD	0.6	Exp. Date	1402/12/29	Lot. No	
High Control	R&D	Mean	34.3	SD	0.8	Exp. Date	1402/12/29	Lot. No	
Low Control	R&D	Mean	25.34	SD	0.32	Exp. Date	1402/12/29	Lot. No	

